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Multimode Faserschalter®

instruction manual

multimode fiberswitch®

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Read carefully before switching on the power! Please also see instructions for safety when using piezoelectric actuators and power supplies!

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**Multimode Fiberswitch
control interface: RS232, USB and BCD
or RS232, BCD and Ethernet**

- | | | | |
|----------------|--------------------------|---------------|--------------------------|
| F-109- | <input type="checkbox"/> | | |
| F-162- | <input type="checkbox"/> | F-142- | <input type="checkbox"/> |
| F-163- | <input type="checkbox"/> | F-143- | <input type="checkbox"/> |
| F-164- | <input type="checkbox"/> | F-144- | <input type="checkbox"/> |
| F-166- | <input type="checkbox"/> | F-146- | <input type="checkbox"/> |
| F-169- | <input type="checkbox"/> | F-149- | <input type="checkbox"/> |
| F-1612- | <input type="checkbox"/> | | |

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Description: multimode fiber switcher

to which this declaration relates is in conformity with the following standards or normative documents:

EN 50082, part 2
EN 55011, Class B
EN 60204-1
EN 61010

The declaration is world-wide valid as the manufacturer's declaration of compliance with the requirements of the above mentioned national and international standards.

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1. introduction

This manual describes the fiber switch® from **piezosystem jena**. You will also find additional information regarding piezoelectric products.

If you have any problems please contact the manufacturer of the system: **piezosystem jena GmbH**, Stockholmer Strasse 12, 07747 Jena.
Phone: +49 36 41 66 88-0

2. certification of **piezosystem jena**

The company, **piezosystem jena GmbH**, has worked according to a DIN EN ISO 9001 certified quality management system since 1999. Its effectiveness is verified and proven by periodic audits by the TÜV.

This instruction manual includes important information for using piezo actuators. Please take the time to read this information. Piezo positioning systems are mechanical systems that are capable of the highest precision. Correct handling guarantees the maintenance of this precision over a long period of time.

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Jena, 17. November 2010

4. safety instructions

Please check the completeness of the delivery after receiving the shipment:

- fiber switch
- USB-, RS232- or BCD-cable depending on the fiber switch type
- instruction manual
- CD-ROM with drivers

5. safety instruction

- Do not open the units! The warranty will be lost when the unit is opened!
- There are no user serviceable parts inside and opening or removing covers may expose you to dangerous shock hazards or other risks. Refer all servicing to qualified service personnel.
- Do not spill any liquids into the connectors or use the units near water.
- Do not place any heavy objects on any cables (e.g. power cords, optical cables). Damage may cause malfunction, shock, or fire! Do not bend the optical cables.
- Do not place the units on a sloping or unstable cart, stand, or table as they may fall or not work accurately.
- Only work with the units in a clean and dry environment! Do not use the switch under rough environmental conditions; they are designed for lab applications only.
- Please only use original parts from **piezosystem jena**.
- **piezosystem jena** does not give any warranty for damages or malfunction caused by additional parts not supplied by **piezosystem jena**. Additional cables or connectors will change the specified data. This can change the specified properties of the units and cause them to malfunction.
- Fiberswitches are sensitive systems. They will demonstrate their excellent properties only if they are handled correctly! Please mount them properly, only at the special mounting points.
- Please contact **piezosystem jena** or your local representative, if there are any problems with your fiberswitch.

Immediately unplug your unit from the power supply and refer servicing to qualified service personnel under the following conditions:

- when the power supply cord, connector, or optical cable is damaged
- if the unit has been exposed to rain or water
- if the unit has been dropped or the housing is damaged

6. maintenance and inspection

Before cleaning the exterior box of the voltage amplifier, turn off the power switch and unplug the power plug. Failure to do so may result in a fire or electrical shock.

Clean the exterior box using a damp cloth that has been firmly wrung-out. Do not use alcohols, benzene, paint thinner, or other inflammable substances. If flammable substances come into contact with an electrical component inside the voltage amplifier, it may result in a fire or electrical shock.

7. instructions for checking

When you open the package, please check to make sure all the necessary parts are complete (see packing list) and nothing is damaged.

If there is any damage to the system please contact our local representative immediately! If the packaging material is damaged please confirm this with the shipping company.

8. overview

Depending on the type of the fiber, the switches are suitable for the spectral range of 180 nm...2600 nm. The insertion losses are in the range of 1.4...3.5 dB.

The fibers are moved utilizing piezo-bending actuators. Typical switching time is 2...10 ms.

A complete switch consists of casing, fibers, and fiber connectors. **piezosystem jena** is able to deliver switches equipped with the following standard fiber types:

- 100/140 μm
- 400/440 μm
- 600/720 μm

The control voltage must be in the range determined in the technical data. Any other values will cause the switch to work improperly.

9. technical data of fiber switches

9.1. casing

fiber core/cladding diameter	≤400/440 μm 600/720 μm	100/140μm 400/440μm 600/720μm	400/440μm 600/720μm
fiberswitch	1x2,1x3	1x4,1x6,1x9	1x12
length (mm)	375,5	375.5	375.5
height (mm)	88.9	88.9	88.9
width (mm)	448	448	448

table 1: dimensions of the casing

9.2. electrical parameters

Interface depending on fiber switch type	serial RS232 (SubD 9 pol), USB or Ethernet and BCD
power	100-240 V, fuse 2 A slow
power input	2 W
switching time	2...10 ms
display	one LED per channel

table 2: electrical parameters

9.3. optical parameters

	λ: 400 nm	λ: 800 nm
max. beam power	0.05 W	2.00 W
connectors	SMA – style fiber optic connector (bulkhead)	
adapter fiber**	1 m ±10% connector SMA-SMA*	

table3: optical parameters

Please note: You will find additional optical parameters on the measurement protocol.

9.4. environmental conditions

operating temperature	0°C...+60°C
storage temperature	-40°C...+80°C
humidity	max. 55%

table4: environmental conditions

* other fiber lengths on request

** LWL adapter on request

10. set up procedure with RS232, USB, or Ethernet

RS232:

Please set the parameters for the RS232 as following:

- 115200 baud
- software-handshake
- no parity
- 8 databits
- 1 stop bit

Please use a normal serial cable (1:1, plug/socket) to connect your computer to the RS232 on the fiber switch. Use the installation software to install the driver.

USB:

Parameters for the USB interface are not required, a USB cable is included. An installation program to set up the USB driver is located on the CD or on our website under www.piezosystem.com "products & services."

Ethernet:

With the help of USB - Ethernet converter SX-3000GB the corresponding channels are selected via Ethernet interface. In this option USB is not available. For software installation and instructions, use the included CD.

10.1. interface control:

After powerup the RS232 sends a welcome string to the pc:

"piezosystem jena FSM12" followed by<LF><CR>

"piezosystem jena FSM13" followed by<LF><CR>

"piezosystem jena FSM14" followed by<LF><CR>

"piezosystem jena FSM16" followed by<LF><CR>

"piezosystem jena FSM19" followed by<LF><CR>

"piezosystem jena FSM112" followed by<LF><CR>

10.2. control commands:

<CR> starts command recognition, erases input buffer

ch,1<CR> binary output 1

ch,2<CR> binary output 2

ch,3<CR> binary output 3

ch,4<CR> binary output 4

and so on

11. BCD code

11.1. initiation

The setup procedure is for F-162, 163, 164, 166, 169, 1612, F-142, 143, 144, 146, and 149.

The switching procedure will be realized by BCD code.

The parallel port (LPT1/LPT2) of your computer must be connected with 25 pin D-SUB port of the switch.

If the red LED is illuminated, the switch is in operation mode.

The Software is provided on the CD or on the internet www.piezosystem.com.

11.2. pinning and BCD code table

The TTL-Interface is realized via the D-SUB25pin-socket. The cable is rated at a 2 m cable length.

D-SUB 25pin-socket	function
2	D0
3	D1
4	D2
5	D3
12, 13	+5 V / 250 mA out
24, 25	GND

table 5: pinning

D3	D2	D1	D0	fiber switch					
				1x2	1x3	1x4	1x6	1x9	1x12
0	0	0	0	1	1	1	1	1	1
0	0	0	1	2	2	2	2	2	2
0	0	1	0		3	3	3	3	3
0	0	1	1			4	4	4	4
0	1	0	0				5	5	5
0	1	0	1				6	6	6
0	1	1	0					7	7
0	1	1	1					8	8
1	0	0	0					9	9
1	0	0	1						10
1	0	1	0						11
1	0	1	1						12

table 6: BCD-code

12. your notes