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10. Ihre Notizen

table of contents

1.	introduction.....	17
2.	certification of <i>piezosystem jena</i>	17
3.	Declaration of conformity	17
4.	purchased part package	18
5.	instructions for using piezoelectric elements and power supplies.....	18
6.	safety instructions.....	19
6.1	maintenance and inspection	21
6.2	environmental conditions	21
7.	instructions for checking the function of the system / quick start	21
8.	description of the voltage amplifier: 30V300	22
8.1.	common introduction	22
8.2.	technical data	23
8.3.	initiation.....	24
8.4.	service.....	24
8.5.	modulation input: MOD	25
8.6.	monitor output socket: MON	25
8.6.1.	30V300 / 30V300nanoX:.....	25
8.6.2.	30V300CLE / 30V300nanoX CLE:	25
8.7.	actuator-socket: OUT.....	25
8.8.	closed-loop button (CLE only).....	26
8.9.	ERROR-LED	26
9.	possibilities of the error correction.....	26
10.	your notes	27

1. introduction

This manual describes the 30V300 voltage amplifier series from **piezosystem jena**. You will also find additional information regarding piezoelectric products.

Definition: All systems from **piezosystem jena**, such as electronics, actuators, and optical systems, are called units.

If you have any problems please contact the manufacturer of the system: **piezosystem jena GmbH**, Stockholmer Strasse 12, 07747 Jena.

Phone: +49 36 41 66 88-0

2. certification of *piezosystem jena*

The company **piezosystem jena GmbH** has worked according to a DIN EN ISO 9001 certified quality management system since 1999. Its effectiveness is verified and proven by periodic audits by the TÜV.

This instruction manual includes important information for using piezo actuators. Please take the time to read this information. Piezo positioning systems are mechanical systems that are capable of the highest precision. Correct handling guarantees the maintenance of this precision over a long period of time.

3. Declaration of conformity

The CE Declaration of Conformity is available on our website: www.piezosystem.com.

4. purchased part package

Please check the completeness of the delivery after receiving the shipment:

- 30V300 piezo amplifier
- wide range power supply 24 V DC
- instruction manual
- CD-ROM with driver, software, and instruction manual

5. instructions for using piezoelectric elements and power supplies

- Piezoelectric actuators from **piezosystem jena** are controlled by voltages of up to 150 V. These values can be quite hazardous. Therefore, read the installation instructions carefully and ensure that only authorized personnel handle the power supply.
- After transportation, piezoelectric actuators should be allowed to adapt to room temperature for approximately 2 hours before being switched on.
- Piezoelectric actuators are made from ceramic materials, with and without metallic casings. The piezo-ceramic is a relatively brittle material. This should be noted when handling piezoelectric actuators. All piezo elements are sensitive to bending or shock forces.
- Due to the piezoelectric effect, piezoelectric actuators can generate electrical charges by changing the mechanical load or the temperature, or by such actions as those described above.
- Piezoelectric actuators are able to work under high compressive forces, but only actuators with pre-load can be used under tensile loads (these tensile forces must be less than the pre-load, given in the data sheet).
- Please note that the acceleration of the ceramic material (e.g., caused by fall down, discharging, or high dynamic application) can cause damage to the actuator.
- Heating of the ceramic material will occur during dynamic operation and is caused by structure conditional loss processes. This may cause failure if the temperature exceeds the specified values cited below.
- When temperature increases up to the Curie temperature (usual values approx. 140°C - 250°C), the piezoelectric effect disappears.

- Piezoelectric actuators, such as stacks or various tables, work electrically as a capacitor. These elements are able to store electrical energy over a long period of time, (up to some days) and the stored energy may be dangerous.
- If the actuator remains connected to the drive electronics, it is discharged within a second after shutdown and quickly reaches harmless voltage values.
- Piezoelectric actuators can only generate voltages by warming or cooling. The discharge potential should not be ignored due to the inner capacitance. This effect is insignificant at room temperature.
- Piezo actuators from **piezosystem jena** are adjusted and glued. Any opening of the unit will cause misalignment or possible malfunction, and will make the guarantee invalid.
- Please contact **piezosystem jena**, or your local representative, if there are any problems with your actuator or power supply.

Caution! Shock forces may damage the built-in ceramic element. Please avoid such forces, and handle the units with care, otherwise the guarantee will become invalid.

Please note that after a long period of storage, the piezo stack in your piezo actuator might experience a certain extent of depolarization. Your amplifier is equipped with a soft start function, which, after switching on, slowly goes through the whole motion range once. This corrects any depolarization effects which might have occurred.

- If your application does not allow such a procedure, please let us know.
- This function can only be deactivated at **piezosystem jena**.

6. safety instructions

Icons:

RISK OF ELECTRIC SHOCK! Indicates that a risk of electric shock is present and the associated warning should be observed.

CAUTION! REFER TO OPERATOR'S MANUAL – Refer to your operator's manual for additional information, such as important operating and maintenance instructions.

RISK OF ELECTRIC SHOCK!

- Do not open the units! There are no user serviceable parts inside, and opening or removing covers may expose you to dangerous shock hazards or other risks. Refer all servicing to qualified service personnel.
- Do not spill any liquids into the cabinet or use the units near water.

CAUTION!

- Allow adequate ventilation around the units so that heat can properly dissipate. Do not block ventilated openings or place the units near a radiator, oven, or other heat sources. Do not put anything on top of the units except those that are designed for that purpose (e.g. actuators).
- Only work with the units in a clean and dry environment! Only specially prepared units (e.g. actuators) can work under other conditions!
- Please only use original parts from **piezosystem jena**. **piezosystem jena** does not give any warranty for damages or malfunction caused by additional parts not supplied by **piezosystem jena**. Additional cables or connectors will change the calibration and other specified data. This can change the specified properties of the units and cause them to malfunction.
- Piezo elements are sensitive systems capable of the highest positioning accuracy. They will only demonstrate their excellent properties if they are handled correctly! Please mount them properly only at the special mounting points.
- Do not insert objects of any kind into the cabinet slots, as they may touch dangerous voltage points, which can be harmful or fatal, or may cause electric shock, fire, or equipment failure.
- Do not place heavy objects on any cables (e.g. power cords, sensor cables, actuator cables, optical cables). Damage may cause malfunction, shock, or fire!
- Do not place the units on a sloping or unstable cart, stand, or table as they may fall or not work accurately.

Immediately unplug your unit from the wall outlet and refer servicing to qualified service personnel under the following conditions:

- when the power supply cord or plug is damaged
- if liquid has been spilled or objects have fallen into the unit
- if the unit has been exposed to rain or water
- if the unit has been dropped or the housing is damaged

6.1 maintenance and inspection

CAUTION!

- Before cleaning the exterior box of the voltage amplifier, turn off the power switch and unplug the power plug. Failure to do so may result in a fire or electrical shock.
- Clean the exterior box using a damp cloth that has been firmly wrung-out. Do not use alcohols, benzene, paint thinner, or other inflammable substances. If flammable substances come into contact with an electrical component inside the voltage amplifier, it may result in a fire or electrical shock.

6.2 environmental conditions

The amplifier can be used:

- indoors only
- at an altitude up to 2000 m
- at a temperature of 5 ... 35 °C
- at a relative humidity of 5 ... 95% (non-condensing)

The recommended environmental conditions:

- indoors only
- at an altitude up to 2000 m
- at a temperature of 20...22 °C
- at a relative humidity of 5 ... 80% (non-condensing)

7. instructions for checking the function of the system / quick start

When you open the package, please check to make sure all the necessary parts are complete (see packing list) and nothing is damaged.

Check the electronics and the actuator for any visible damage:

- The top and bottom plates of the actuator (if it does not have another shape) should be parallel each to each other, and not contain any scratches.
- If there is any damage to the system please contact our local representative immediately!
- If the packaging material is damaged please confirm this with the shipping company.

quick start:

- The power switch is OFF.
- Check that the available main voltage is the same as the label on the amplifier's rear panel.
- Connect the power supply.
- Connect the piezo element by using the "OUT" connector.
- Connect the measurement system on the "SENSOR" connector, if available.
- Switch on the system by using the POWER switch on the front panel.
- The green "ON" LED is permanently lit. A short self test of all parts will follow.
- Whenever the system is switched on, it slowly goes through the voltage range from -20 V to 130 V. This is to protect the actuator against depolarization effects.
- Switch on the closed loop by pressing the closed loop button. The "CL" LED will light up.
- Now you can set the position with the OFFSET potentiometer.
- If the "OOR" LED lights up, there is an error.
- Please switch off the device.
- The function check is complete.

8. description of the voltage amplifier: 30V300

8.1. common introduction

This amplifier is designed for use as a single unit in industrial settings. It is compact, robust, and mountable in different ways, as well as highly reliable. The 30V300 series was designed for universal use with a wide main supply voltage from 10 V to 30 V DC. The casing is available as a screw slot version or for mounting in a 19" rack mount casing.

Piezo actuators can be controlled in two different ways: by using the potentiometer on the front panel or by an analog signal between 0 and +10 V.

A unique feature of the 30V300 is that it can be used in combination with strain gauge or capacitive feedback sensors without additional modification.

The 30V300 amplifier contains a temperature surveillance circuit. Whenever the internal temperature rises above a preset limit, the active cooling is activated. If, due to an abnormality or malfunction, the temperature increases even more, the

amplifier output is disabled until the temperature decreases again. After cooling down the amplifier will work normally.

All operating elements are comfortably operated from the front side. Extension cables for great distances are available at **piezosystem jena**.

8.2. technical data

supply voltage [V]	10V... 30 DC
max. current	5A @ 10V DC (2A @ 24V DC)
power supply socket	socket with output connector 2,1mm
dimensions w x h x d [mm]	130mm x 86mm x 230mm
mass [g]	1800
channels	1
output power [W]	39
output current [mA]	300 (constant) 2 * 150 (nanoX only)
output voltage [V]	-20 ... 130
output voltage 2 [V]	130 ... -20 [nanoX only]
output connector (OUT)	LEMO 0S.302 (30V300 / 30V300CLE) ODU 3pin (30V300nanoX / 30V300nanox CLE)
modulation input (MOD)	0...10V BNC
monitor voltage (MON)	-2 ... +13V BNC (0...10V using CLE)
inner resistance	<100kΩ OL / <1kΩ CLE
output noise	<0,3 mV _{RMS} @ 500Hz
bandwidth [Hz]	4000
polarity	positiv
measuring system connector (CLE only)	ODU 4pin
features	temperature controlled fan, short circuit proof

table 1: technical data

8.3. initiation

Please connect the device with the wall outlet by using the external power supply. Please make sure that it is suitable to the type of amplifier being used. The LED ON lights up after the power is switched on. If using the modulation signal and monitor signal, pay attention to the correct assignment of them on the front panel BNC sockets. The actuators are driven by voltages of up to +130 V. Please pay attention to shock hazard protection.

After power up, the soft start function is started for 10 seconds to initiate the actuator. During this routine, the actuator is moved once through the whole voltage or motion range and back. While in soft start, the modulation input and the offset potentiometer are disabled. After the soft start, the amplifier will be working normally and the PID controller will be working in the state which was active before power off.

A special protection circuit is able to detect shorts of the actuator or the cables. In case of a short circuit, the amplifier stages are disabled. Once every 4 seconds the device checks the output for 200 ms to determine whether or not the short is still present. If the short is removed the amplifier continues to work normally.

8.4. service

The actuator is connected to the amplifier by a LEMO socket (ODU 3pin for nanoX type). The potentiometer is used to adjust the level of the output voltage (OFFSET) in open loop mode. In closed loop mode, the actuator position can be adjusted using the potentiometer. The rest position of the actuator is hereby adjusted. An external triggering signal in the range of 0 to +10 V (MOD input socket) can control the output voltage. It is possible to use this input for highly dynamic scanning functions. To avoid damage to the actuators, it is recommended that the potentiometer be adjusted to its maximum counter clockwise position before switching on the amplifier.

By superposition of the modulating voltage and adjusted offset, voltages of up to 150 V might occur. This operating state should be avoided to protect the ceramics and increase the MTBF. In addition, the OOR-LED lights up in the case of a limit error and there is no automatic disconnection of the overvoltage when this occurs!

We recommend switching the amplifier on approx. 2 hours before measurement in the sub- μm range takes place to guarantee stable temperature circumstances. A constant temperature environment is necessary for precise positioning tasks. Please note that a temperature change of $\Delta T = 5 \text{ K}$ will cause a $13 \mu\text{m}$ increase in length of a 20 cm steel rod.

Special qualities of piezo ceramics, like hysteresis and creep, can cause inaccuracies in the case of non-observance. These basic qualities of piezo elements are described in the "piezoline" tutorial in our catalogue. Do not hesitate to contact our staff if you need more information.

After switching ON, the green LED indicates "ready to operate". If the OFFSET potentiometer is in the maximum counter clockwise position, the actuator is in the rest position. By turning the potentiometer into the maximum clockwise position, the maximum motion will be done. The motion depends on the actuator specification. Please refer the hysteresis data curve for details.

The red OOR LED indicates overload or underload. Please avoid over- and underload in these cases or switch off the device, if the LED still flashes.

8.5. modulation input: MOD

The motion of the actuator may be remotely controlled by using this input. The control signal must be in the range of 0 to +10 V, which can be reached in various ways, such as by using a function generator. There is an internal addition of the MOD signal and the adjusted OFFSET potentiometer.

8.6. monitor output socket: MON

8.6.1. 30V300 / 30V300nanoX:

The output voltage in ratio 10:1 (-2 to +13 V) is available at this socket and can be monitored by using an analogue oscilloscope. This is worth doing, especially during dynamic operation.

Pay attention to the inner resistance of the monitor output, only high ohmic measurement devices should be used. The output is short protected, please avoid injected voltages.

8.6.2. 30V300CLE / 30V300nanoX CLE:

The measuring voltage (0...10 V) is available at this socket and can be monitored by using an analogue oscilloscope. This is worth doing, especially during dynamic operation.

Pay attention to the inner resistance of the monitor output, only high ohmic measurement devices should be used. The output is short protected, please avoid injected voltages.

8.7. actuator-socket: OUT

Please connect the actuator by using this socket.

8.8. closed-loop button (CLE only)

The 30V300CLE/ 30V300nanoX CLE amplifiers include an electronic closed-loop system. This compensates for the drift and hysteresis of the piezo actuator. Please allow for settling time when using in dynamic operation.

8.9. ERROR-LED

The ERROR-LED signals failure of an operating voltage or it signals that an operating voltage is not in the acceptable range. Possible causes can be over modulation or an inadequate external current supply.

9. possibilities of the error correction

Please protect the power cords from harm or shorts, which could cause malfunction of the equipment. Large strains on the cords can cause interruptions. When using the open loop configuration, the output voltage can be read from the display without plugging in the actuator.

error	possible correction
ERROR-LED is illuminated	Please check the power supply.
erroneous, illogical display values	Please check the sensor cable and connection.
ON-LED doesn't work	Please check the power supply and the power cord.
OOR-LED is illuminated	Please check the position of the adjustable potentiometer and the amplitude of the modulation signal. High modulating voltage overdrives the amplifier. Decrease the modulation signal or reduce the control voltage until OOR-LED is off.

table2: possible error correction

The equipment is calibrated before delivery for the appropriate actuator with integrated measuring system. An exchange of the actuator-amplifier-combination leads to inaccuracies in setting movement and loss of positioning accuracy. If the equipment is damaged during the use of another actuator, no assurance can be undertaken!

Equipment customization makes adaptations, according to customer preferences, that are possible in terms of the technical threshold values e.g. the main voltage or the output voltage. Please contact our technical service department in order to find out the possibilities for your specific application. Special adaptations must be paid for by the customer.

10. your notes