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9. Ihre Notizen

Table of Contents

| | | |
|------|--|----|
| 1. | Introduction | 21 |
| 2. | Certification of piezosystem jena | 21 |
| 3. | Declaration of conformity | 22 |
| 4. | Instructions for using piezoelectrical elements and power supplies | 23 |
| 5. | Safety instructions | 24 |
| 5.1. | maintenance and inspection | 25 |
| 5.2. | environmental condition | 25 |
| 6. | Instructions for checking the function of the system / quick start | 26 |
| 7. | How to operate the 12V40 | 27 |
| 7.1. | Technical data | 27 |
| 7.2. | Initialization | 28 |
| 7.3. | Service | 28 |
| 7.4. | SMB modulation input: MOD | 29 |
| 7.5. | SMB monitor output socket: MON | 29 |
| 7.6. | Actuator socket: OUT | 29 |
| 7.7. | Switch "closed loop" (12V40SG / 12V40CLE) | 29 |
| 7.8. | Wiring of the rear panel connector | 30 |
| 8. | Possibilities of the error correction | 32 |
| 9. | Your notes | 33 |

1. Introduction

This manual describes the 12V40 voltage amplifier series from **piezosystem jena**. This series contains the following amplifiers: 12V40, 12V40SG, 12V40CLE, 24V40, 24V40SG, 24V40CLE, 12V40C and 12V40CSG. This manual also offers you additional information about piezoelectric products.

Definition: All systems from **piezosystem jena**, such as electronics, actuators, and optical systems, are called units.

If you have any problems please contact the manufacturer of the system: **piezosystem jena GmbH**, Stockholmer Strasse 12, 07747 Jena.
Phone: +49 36 41 66 88-0

2. Certification of **piezosystem jena**

The company **piezosystem jena GmbH** has worked according to a DIN EN ISO 9001 certified quality management system since 1999. Its effectiveness is verified and proven by periodic audits by the TÜV.

The latest information on our products and corresponding certificates can be found on our website <http://www.piezosystem.de> or requested from our sales team.

This instruction manual includes important information for using piezo actuators. Please take the time to read this information. Piezo positioning systems are mechanical systems that are capable of the highest precision. Correct handling guarantees the maintenance of this precision over a long period of time.

3. Declaration of conformity

You can request the declaration of conformity from our sales team.

4. Instructions for using piezoelectrical elements and power supplies

- Piezoelectric actuators from **piezosystem jena** are controlled by voltages of up to 150 V. These values can be quite hazardous. Therefore, read the installation instructions carefully. Only authorized personnel should handle the power supply.
- After transportation, piezoelectric actuators should be allowed to adapt to room temperature for approximately 2 hours before being switched on.
- Piezoelectric actuators are made from ceramic materials with and without metallic casings. The piezo-ceramic is a relatively brittle material. This should be noted when handling piezoelectrical actuators. All piezo elements are sensitive to bending or shock forces.
- Due to the piezoelectric effect, piezo actuators can generate electrical charges by changing the mechanical load or the temperature, or by such actions as those described above.
- Piezoelectric actuators are able to work under high compressive forces, but only actuators with pre-load can be used under tensile loads (these tensile forces must be less than the pre-load, given in the data sheet).
- Please note that the acceleration of the ceramic material (e.g., caused by fall down, discharging, or high dynamic application) can cause damage to the actuator.
- Heating of the ceramic material will occur during dynamic operation and is caused by structure conditional loss processes. This may cause failure if the temperature exceeds the specified values cited below.
- With increasing temperature, up to the Curie temperature (usual values approx. 140°C - 250°C), the piezoelectric effect disappears.
- Piezoelectric actuators, such as stacks or various tables, work electrically as a capacitor. These elements are able to store electrical energy over a long period (up to some days) and the stored energy may be dangerous.
- If the actuator remains connected to the drive electronics, it is discharged within a second after shutdown and quickly reaches harmless voltage values.
- Piezo actuators can only generate voltages by warming or cooling. The discharge potential should not be ignored due to the inner capacitance. This effect is insignificant at usual room temperature.
- Piezo actuators from **piezosystem jena** are adjusted and glued. Any opening of the unit will cause misalignment or possible malfunction, and the guarantee will become invalid.
- Please contact **piezosystem jena**, or your local representative, if there are any problems with your actuator or power supply.

Caution! Shock forces may damage the built-in ceramic element. Please avoid such forces, and handle the units with care, otherwise the guarantee will become invalid.

5. Safety instructions

Icons:

RISK OF ELECTRIC SHOCK! Indicates that a risk of electric shock is present and the associated warning should be observed.

CAUTION! REFER TO OPERATOR'S MANUAL – Refer to your operator's manual for additional information, such as important operating and maintenance instructions.

RISK OF ELECTRIC SHOCK!

- Do not open the units! There are no user serviceable parts inside. Opening or removing covers may expose you to dangerous shock hazards or other risks. For servicing please refer to qualified service personnel.
- Do not spill any liquids into the cabinet or use the units near water.

CAUTION!

- Allow adequate ventilation around the units so that heat can properly dissipate. Do not block ventilated openings or place the units near a radiator, oven, or other heat sources. Do not put anything on top of the units except accessories designed for that purpose (e.g. actuators).
- Only work with the units in a clean and dry environment! Only specially prepared units (e.g. actuators) can work under other conditions!
- Please only use original parts from **piezosystem jena**. **piezosystem jena** does not give any warranty for damages or malfunction caused by additional parts not supplied by **piezosystem jena**. Additional cables or connectors will change the calibration and other specified data. This can change the specified properties of the units and cause them to malfunction.
- Piezoelements are sensitive systems capable of the highest positioning accuracy. They will only demonstrate their excellent properties if they are handled correctly! Please mount them properly, only at the special mounting points.
- Do not insert objects of any kind into the cabinet slots, as they may touch dangerous voltage points, which can be harmful or fatal or may cause electric shock, fire, or equipment failure.
- Do not place heavy objects on any cables (e.g. power cords, sensor cables, actuator cables, optical cables). Damage may cause malfunction, shock, or fire!
- Do not place the units on a sloping or unstable cart, stand, or table as they may fall or not work accurately.

Immediately unplug your unit from the wall outlet and refer servicing to qualified service personnel under the following conditions:

- when the power supply cord or plug is damaged
- if liquid has been spilled or objects have fallen into the unit
- if the unit has been exposed to rain or water
- if the unit has been dropped or the housing is damaged

5.1. Maintenance and inspection

CAUTION!

- Before cleaning the exterior box of the voltage amplifier, turn off the power switch and unplug the power plug. Failure to do so may result in a fire or electrical shock.
- Clean the exterior box using a damp cloth that has been firmly wrung-out. Do not use alcohols, benzene, paint thinner, or other inflammable substances. If flammable substances come into contact with an electrical component inside the voltage amplifier, it may result in a fire or electrical shock.

5.2. Environmental condition

The amplifier can be used:

- indoors only
- at an altitude up to 2000 m
- at a temperature of 5 ... 35 °C
- at a relative humidity of 5 ... 95% (non-condensing)

The recommended environmental conditions:

- indoors only
- at an altitude up to 2000 m
- at a temperature of 20...22 °C
at a relative humidity of 5 ... 80% (non-condensing)

6. Instructions for checking the function of the system / quick start

When you open the package, please check to make sure all the necessary parts are complete (see packing list) and nothing is damaged.

Check the electronics and the actuator for any visible damage:

- The top and bottom plates of the actuator (if it does not have another shape) should be parallel each to each other, and not contain any scratches.
- If there is any damage to the system, please contact our local representative immediately!
- If the packaging material is damaged please confirm this with the shipping company.

Before you switch on the system, please check:

- that the main voltage supplied in your country is the same as the one installed for the external power supply (not included in the shipment!).
- that the secondary voltage (DC) matches with the amplifier.
- Connect the power cable.
- Connect the piezo element by using the LEMO connector.
- Turn off closed loop (12V40SG / 12V40CLE).
- Switch on the power supply. The red LED indicates that the amplifier is working.
- By turning the potentiometer to its maximum position, the actuator will reach its maximum movement.
- Offset is accomplished by the offset potentiometer or external offset (e.g. a load). Because of the fast motion of the actuator a noise or crack can be heard. This is normal and not a malfunction.
- Turn the potentiometer into its minimum position.

12V40SG / 12V40CLE:

- Turn on closed loop. Turn the potentiometer into the maximum position so that the actuator does its maximum closed loop movement. Turn the potentiometer into its minimum position and switch off closed loop.
- Switch off the power supply. Ensure that the potentiometer is in its minimum position.
- The function check is now complete.

7. How to operate the 12V40

The 12V40 voltage amplifier series was especially developed for one channel positioning tasks in optics, laser physics, microbiology, and machining. With an output noise less than 0.3 mV, it is well suited for positioning in the sub-nm range. It is also well suited for OEM applications.

7.1. Technical data

| | |
|----------------------------|---|
| supply voltage | +12V DC (optional +24V DC) |
| max. current | 1A @ 12V DC (0.5A @ 24V DC) |
| dimensions BxHxT [mm] | 173mm x 130mm/105mm x 45mm |
| mass [g] | 800 |
| channels | 1 |
| output power [W] | 6 |
| output current [mA] | 40 |
| output voltage (OUT) | -10...+150 V |
| output connector (OUT) | LEMO 0S.302 |
| modulation input (MOD) | 0...+10V SMB |
| inner resistance | 10k Ω |
| monitor voltage (MON) | -2...+15V SMB -2...+15V rear panel |
| Inner resistance | 100k Ω (SMB) / 10k Ω (rear panel) |
| output noise | <0,3 mV _{RMS} @ 500Hz |
| bandwidth [Hz] | 400 |
| polarity | positiv |
| measuring system | strain gauges (12V40SG) external (12V40CLE) |
| measuring system connector | LEMO 0S.304 (12V40SG) ODU 4pin series L (12V40CLE) |

table 1: technical data 12V40 series

7.2. Initialization

Please connect the device to the wall outlet by using the external power supply. It must be suitable for the type of amplifier being used. The power LED ON lights up after the power supply is switched on. If using the modulation signal and monitor signal, make sure that the correct front panel SMB sockets are used.

The actuators are driven by voltages of up to +150 V. Please pay attention to shock hazard protection.

7.3. Service

The actuator is connected to the amplifier by a LEMO socket. The potentiometer is used to adjust the level of the output voltage (OFFSET). The rest position of the actuator is hereby adjusted. An external triggering signal in the range of 0...+10 V (MOD input socket) can control the output voltage. It is possible to use this input for highly dynamic scanning functions. To avoid damage to the actuators, it is recommended that the potentiometer be adjusted to its maximum counter clockwise position before switching on the amplifier.

By superposition of the modulating voltage and adjusted offset, voltages up to 160 V might occur. This operating state is to be avoided to protect the ceramics and increase the MTBF. In addition, the OVL-LED lights up in the case of a limit error and there is no automatic disconnection of the overvoltage when this occurs! To guarantee stable temperature circumstances, we recommend that the amplifier be turned on approx. 2 hours before measurement in the sub-nm range takes place. A constant temperature environment is necessary for precise positioning tasks. Please note that a temperature change of $\Delta T = 5 \text{ K}$ will cause a 13 μm increase in 20 cm length of steel rod.

The special qualities of piezo-ceramics, such as hysteresis and creep, can cause inaccuracies in the case of nonobservance. These basic qualities of piezo-elements are described in the "piezoline" tutorial in our catalog. Do not hesitate to contact our staff if you need further information.

After switching ON, the red LED indicates "ready to operate". If the OFFSET potentiometer is in its maximum counter clockwise position, the actuator is in its rest position. By turning the potentiometer into the maximum clockwise position, the maximum motion will be done. The motion depends on the actuator's specifications so please refer to the hysteresis data curve for details.

The red OVL-LED indicates overload, while the red UDL-LED indicates underload. Please avoid over- & underload in these cases or switch off the device if the LED is still lit.

7.4. SMB modulation input: MOD

The motion of the actuator may be remotely controlled by using this input. The control signal must be in the range of 0 ... +10 V. There is an internal addition of the MOD signal and the adjusted OFFSET potentiometer. If there is an HIGH TTL level at the Pin8 of the rear panel connector, OFFSET and MOD (modulation input) are switched off. The remote control via Pin7 of the rear panel connector is active and possible.

7.5. SMB monitor output socket: MON

12V40:

The output voltage in ratio 10:1 (-1 ... +15 V) is available at this socket and can be monitored by using an analog oscilloscope. This is especially recommended for dynamic operation. Pay attention to the inner resistance of the monitor output since only high ohmic measurement devices should be used. The output is short protected so please avoid injected voltages.

12V40SG / 12V40CLE:

The voltage of the position sensor is displayed via this socket. The value is 0 to 10V nominal in closed-loop operation.

In case of mechanical overshoot of the used actuator and in open-loop operation, the voltage can assume values from -15 to 15V.

The output signal can be monitored by using an analog oscilloscope. This is especially recommended for dynamic operation. Pay attention to the inner resistance of the monitor output since only high ohmic measurement devices should be used. The output is short protected so please avoid injected voltages.

7.6. Actuator socket: OUT

Please connect the actuator by using this socket.

7.7. Switch “closed loop” (12V40SG / 12V40CLE)


Closed loop operation can be turned on by using this switch. The 12V40SG amplifier includes an electronic closed loop and a strain gauge measurement system. The 12V40CLE uses a measuring system which is combined with the actuator. This allows for compensation of drift and hysteresis of the piezo actuator. Please take care of settling time when using in dynamic operation.

7.8. Wiring of the rear panel connector

To realize OEM applications, all main control signals are also available at the rear panel connector.

| Pin | name | description |
|-----|------|---|
| 1,4 | GND | ground (Signal) |
| 2 | VCC | +5 V output, for switching SRC |
| 6 | MON | monitor |
| 7 | MOD | modulation input |
| 8 | SRC | source, choice of the amplifier control by HIGH or LOW TTL level HIGH: Pin 7 is active, OFFSET potentiometer and MOD (modulation input) at the front panel are switched off LOW: Pin 7 is not active, motion of the actuator will be controlled by the OFFSET potentiometer and the MOD (modulation input) at the front panel |
| 12 | OUT | output voltage -10...+150 V (* option) |
| 15 | GND | ground (operating voltage) |
| 16 | PWR | + operating voltage |

table 2: pin wiring of the rear panel connector (12V40 / SG / CLE)


* option: is set at the factory, must be specified when ordering

The 12V40C/CSG voltage amplifiers do not have a casing. They should be integrated into customer systems for OEM applications.

| Pin | name | description |
|---------|------|---|
| 1,2,8,9 | GND | ground (Signal) |
| 4 | VCC | +5 V output, for switching SRC |
| 13 | MON | monitor signal |
| 14 | MOD | modulation input |
| 16 | SRC | source, choice of the amplifier control by HIGH or LOW TTL level HIGH: Pin 7 is active, OFFSET potentiometer and MOD (modulation input) at the front panel are switched off LOW: Pin 7 is not active, motion of the actuator will be controlled by the OFFSET potentiometer and the MOD (modulation input) at the front panel |
| 24 | OUT | output voltage -10...+150 V (* option) |
| 30 | GND | ground (operating voltage) |
| 32 | PWR | + operating voltage |

table 3: pin wiring of the rear panel connector (12V40C / 12V40CSG)

* option: is set at the factory, must be specified when ordering

8. Possibilities of the error correction

If the system does not work properly, please check the power cord first.

| error | possible correction |
|----------------------|--|
| ON-LED is off | please check the power supply |
| OVL or UVL-LED flash | over- or underload, please check the modulation input and offset |

table 4: possibilities of the error correction

The equipment is calibrated before delivery for the appropriate actuator with integrated measuring system. An exchange of the actuator-amplifier-combination leads to inaccuracies in setting movement and loss of positioning accuracy. If the equipment is damaged during the use of another actuator, no assurance can be undertaken!

The equipment concept makes adaptations to customer preferences concerning the technical threshold values, such as the main voltage or the output voltage, possible. Please contact our technical service department in order to discover the possibilities for your specific problem.

Any adjustments necessary are to be paid for by the customer.

9. Your notes