

MIPOS LFA

ordering example


The piezo objective focus system MIPOS LFA consists of 3 different parts which need to be ordered separately. Please see the example for a MIPOS unit with 500µm travel range and integrated closed loop feed back sensor enclosed. The MIPOS LFA should be mounted on a microscope with a threading size of W0.8x1/26" at the revolver. The used micro objective shows a threading size of M25x0.75mm. Please see the order example below:

Example for the following configuration:

- MIPOS 500 SG (travel range 500µm with strain gage feed back sensor)
- Microscope revolver threading size: W0.8x1/36" (RMS)
- Micro objective lens threading size: M25x0.75

Configuration:

1. Selection Flex Adapter for W0.8x1/36"
part no.: O-304-11
2. Selection of MIPOS 500 SG
part no.: O-350-01-LFA
3. Selection Lens Flex Adapter M25x0.75
part. no. O-303-99



All part numbers are give on page 2 of the data sheet

or at www.piezosystem.com

Specifications subject to change without notice!